

ICARA[®] Radial Head Prosthesis



Bipolar radial head prosthesis

The **ICARA®** system allows for precise intra-operative adjustment and fine-tuning of the prosthesis preload to ensure optimal tension and stability of the joint.



Distalization of the center of rotation

stable fixation by shifting the torque of the prosthesis to the level of the radial tuberosity




Various stem options

cemented and cementless stems with variable collar height for the prevention of stem subsidence



TiN coating

optimization of wettability and minimization of metal ion release at the articulation surface





Adjustable neck length

avoidance of overstuffing or gapping



Multidirectional self-adjustment

minimization of force peaks on the capitulum of the humerus due to optimal alignment of the head



walk with us

implantcast GmbH

Lüneburger Schanze 26
21614 Buxtehude
Germany

tel.: +49 4161 744-0
fax: +49 4161 744-200

info@implantcast.de
www.implantcast.de

CE 0482

