

# CARPOFIT<sup>®</sup> TMC joint prosthesis



# Trapeziometacarpal joint prosthesis

The **CARPOFIT**® trapeziometacarpal joint system is used to treat rhizarthrosis and rheumatoid arthritis. Outstanding long-term results demonstrate the clinical success of the prosthesis.



## Innovative articulating materials

articulation of ceramic coating (TiNbN) with PEEK-C



## Modular

adaption to patient's anatomy



## Bone-saving

less bone loss in the os trapezium due to small cup diameter



## Long-term results<sup>1,2</sup>

15 year follow-up  
95,6% survival rate



## Pressfit cup

good primary stability




## Rotational stability

hexagonal stem design

<sup>1</sup>2023; Joosten et al.; Ergebnisse von 560 implantierten Daumensattelgelenk-Endoprothesen über einen Zeitraum von 15 Jahren

<sup>2</sup>based on MoM-design



walk with us

**implantcast GmbH**

Lüneburger Schanze 26  
21614 Buxtehude  
Germany

tel.: +49 4161 744-0  
fax: +49 4161 744-200

[info@implantcast.de](mailto:info@implantcast.de)  
[www.implantcast.de](http://www.implantcast.de)

CE 0482

